



Simple 3-Step Testing

1

Step 1: Insert the Test Strip

Turn on the meter by inserting the test strip. No Coding required. The Prodigy Voice® speaks the date and time.


 "Thank you for using Prodigy Voice®. Today's date December 9, 2010. The time 4:04 PM. Please apply blood to the test strip."



2

Step 2: Apply Blood

Follow the audible prompts and apply your blood sample. The meter beeps and then counts down from 6 to 1.


 "Now testing."




3

Step 3: Hear Results

Hear and see your accurate results in 7 seconds.

 "Your blood glucose 6.6 millimoles per liter."



Recalling Test Results

The Prodigy Voice® stores the 450 most recent blood glucose test results with date and time in its memory.

Option 1: Individual Averages

Turn the meter on by pressing the "M" button. Press the "M" button again and the meter will display and audibly report the most recent test result. Use the "Up" and "Down" buttons to navigate.




Option 2: Averages

With the meter off, press and hold the "M" button for five (5) seconds. The meter will display and audibly report the 7, 14, 21, 30, 60 and 90-day averages. Use the "Up" and "Down" buttons to navigate.

Confirmation Window

With the test strip inserted, bring the meter close to your finger or other test site and place beside the blood drop. Blood is automatically drawn into the absorbent channel once the test strip is close enough to the blood drop. DO NOT smear blood onto the test strip or place the blood on top of the test strip.



For Customer Care call: 1.800.243.2636

www.prodigymeter.com

PRODIGY®

Prodigy Diabetes Care, LLC
P.O. Box 481928
Charlotte, North Carolina 28269 USA

©2012 Prodigy Diabetes Care, LLC. Prodigy® and Prodigy Voice® are registered trademarks of Prodigy Diabetes Care. All other trademarks are property of their respective owners. All rights reserved. VQR20 Rev.6 07/12

PRODIGY VOICE®

Blood Glucose Monitoring System



Quick Reference Guide